


Propidium Iodide Stain

F23003



Storage

4 °C in the dark

Product Description

Molecular formula	C ₂₇ H ₃₄ I ₂ N ₄
Molecular weight	668.39 g/mol
Appearance	Clear liquid
Cell permeability	Membrane permeant
Excitation/emission	493/636 nm (in aqueous solution) 533/617 nm (when bound to nucleic acids)

Propidium Iodide Stain is a fluorescent vital dye that binds to nucleic acids. Not being able to permeate intact cell membranes, propidium iodide is taken up by nonviable cells and cells with compromised membranes. Once bound to nucleic acids, its fluorescence increases 20-30 fold and causes the cell to fluoresce red.

Propidium Iodide Stain can be used with Acridine Orange Stain (F23002) to assess cell viability with the automated fluorescence cell counters of the LUNA™ family. Viable nucleated cells will fluoresce green and nonviable nucleated cells will fluoresce red. Due to Förster resonance energy transfer (FRET), the propidium iodide signal absorbs the acridine orange signal in nonviable cells, ensuring no double positive results.

Directions for Use

1. Mix:
 - 1 µL Acridine Orange Stain
 - 1 µL Propidium Iodide Stain
 - 18 µL cell sample
2. Count the sample with a compatible LUNA™.

Disclaimer

This product is for research use only.
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website at www.logosbio.com.
© Logos Biosystems 2018. All rights reserved.



HEADQUARTERS

FL 3
28 Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055
South Korea

Tel: +82 (31) 478-4185

USA

7700 Little River Turnpike STE 207
Annandale, VA 22003
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

EUROPE

1 allée Lavoisier 59650 Villeneuve d'Ascq
France


Tel: +33 (0)3 74 09 44 35

www.logosbio.com

LBSM-RD-PI-PI-001 Rev.2

Propidium Iodide Stain

F23003



Storage

4 °C in the dark

Product Description

Molecular formula	C ₂₇ H ₃₄ I ₂ N ₄
Molecular weight	668.39 g/mol
Appearance	Clear liquid
Cell permeability	Membrane permeant
Excitation/emission	493/636 nm (in aqueous solution) 533/617 nm (when bound to nucleic acids)

Propidium Iodide Stain is a fluorescent vital dye that binds to nucleic acids. Not being able to permeate intact cell membranes, propidium iodide is taken up by nonviable cells and cells with compromised membranes. Once bound to nucleic acids, its fluorescence increases 20-30 fold and causes the cell to fluoresce red.

Propidium Iodide Stain can be used with Acridine Orange Stain (F23002) to assess cell viability with the automated fluorescence cell counters of the LUNA™ family. Viable nucleated cells will fluoresce green and nonviable nucleated cells will fluoresce red. Due to Förster resonance energy transfer (FRET), the propidium iodide signal absorbs the acridine orange signal in nonviable cells, ensuring no double positive results.

Directions for Use

1. Mix:
 - 1 µL Acridine Orange Stain
 - 1 µL Propidium Iodide Stain
 - 18 µL cell sample
2. Count the sample with a compatible LUNA™.

Disclaimer

This product is for research use only.
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website at www.logosbio.com.
© Logos Biosystems 2018. All rights reserved.



HEADQUARTERS

FL 3
28 Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055
South Korea

Tel: +82 (31) 478-4185

USA

7700 Little River Turnpike STE 207
Annandale, VA 22003
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

EUROPE

1 allée Lavoisier 59650 Villeneuve d'Ascq
France

Tel: +33 (0)3 74 09 44 35

www.logosbio.com

LBSM-RD-PI-PI-001 Rev.2